







CIP RECIRCULATOR

ART. NO. 550010

GENERAL

-  The KEOFITT CIP RECIRCULATOR is a device to enable chemical cleaning/disinfection of KEOFITT sampling valves.
-  Compatible with all KEOFITT W9 & Sesame sampling valves with hose piece connections.
-  Cleaning agents not included.
-  **CAUTION!** As with any chemical solutions the care, handling and compatibility with the equipment must be considered before undertaking use of the device for sample valve cleaning. Follow safety measures to avoid chemical contact.

FEATURES

-  Installation: Hose piece
-  Inlet / Outlet: Tube fitting

TECHNICAL DATA

Material (without product contact)

- Bottle (Capacity: 1000 ml) Polypropylene (PP)
- Lid Polypropylene (PP)
- Tube fitting Stainless Steel AISI 316L (1.4404)
O-ring SILICONE, RED (#800830)

Net Weight

- kg/lbs 0.466 kg/1.026 lbs

Spareparts

- NOT APPLICABLE

